

The next generation chart projector with its outstanding design.



## Specifications

Charts	42EA	Dimensions	Main Body : 202(W)×364(D)×172(H) mm
Distance of projection	2.5 ~ 8m		Main Body : 202(W)×364(D)×240(H) mm
Chart magnification	30× (when projected at 5m distance)		Remote Controller : 64(W)×196(D)×20(H) mm
Lamp	PACP-7000 : Halogen lamp (6V-20W)	Weights	PACP-7000 : 4.9 kg
	PACP-7000L : LED lamp		PACP-7000L : 4.7 kg
Masks	35EA (1 open mask, 5 horizontal line masks, 8 vertical line masks, 21 single letter masks)		Remote Controller : 160 g (Battery included)
Filter	1 Red/Green Filter	Basic Accessories	Product Body(PACP-7000 or PACP-7000L), Screen
Chart rotation speed	average 0.21 sec (min. 0.092sec, max.0.324sec)		Remote Controller, Table Stand, Dust Cover, Power Cable
Power supply	AC100-240V, 50/60Hz		T wrench, Fuse, Lamp(PACP-7000), Battery
Power consumption	PACP-7000 : 100VA	Optional accessories	Red/Green glasses, Polarization glasses
	PACP-7000L : 60VA		